

RF Ablation System

Less power,
less time, less pain
...bigger ablation



CE 1936

**RF Generator and
arrays needles
ablation system**

RF ABLATION SYSTEM

-  **Dual mode - temperature and power control**
-  **Large screen for parameter display**
-  **Real time temperature control for enhanced safety**

 Less power,
less pain, same necrosis




 The thinnest rf ablation needle,
making a 4 cm perfectly
spherical necrosis

 A wide range of sizes
and models offers
the best choice for the user



- 15 cm and 20 cm coaxial set models allow the use with the whole extrathin needles range
- Ergonomic and anti-reflux design of the handle increases user friendliness and may be positioned using one hand**
- More echogenic electrode enhances the quality of US image
- Sharp diamond beveled tip facilitates the needle penetration**
- Also available with side port on the handle for fluid injection



ORDERING INFORMATION

EXTRATHIN NEEDLES

CAT. NO.	CANNULA DIAMETER	LENGTH	ARRAYS DIAMETER	NUMBER OF ARRAYS
708172	17 G	15 cm	2 cm	8
708173	17 G	15 cm	3 cm	8
7081735	17 G	15 cm	3,5 cm	8
708174	17 G	15 cm	4 cm	8
708372	17 G	20 cm	2 cm	8
708373	17 G	20 cm	3 cm	8
7083735	17 G	20 cm	3,5 cm	8
708374	17 G	20 cm	4 cm	8

EXTRATHIN NEEDLES WITH COAXIAL INTRODUCER SET

CAT. NO.	CANNULA DIAMETER	LENGTH	ARRAYS DIAMETER	NUMBER OF ARRAYS
708172C	15 G	15 cm / 12,5 cm	2 cm	8
708173C	15 G	15 cm / 12,5 cm	3 cm	8
7081735C	15 G	15 cm / 12,5 cm	3,5 cm	8
708174C	15 G	15 cm / 12,5 cm	4 cm	8
708372C	15 G	20 cm / 17,5 cm	2 cm	8
708373C	15 G	20 cm / 17,5 cm	3 cm	8
7083735C	15 G	20 cm / 17,5 cm	3,5 cm	8
708374C	15 G	20 cm / 17,5 cm	4 cm	8

FUSION 9 NEEDLES WITH SIDE INJECTION PORT:

CAT. NO.	CANNULA DIAMETER	LENGTH	ARRAYS DIAMETER	NUMBER OF ARRAYS
709162	16 G	15 cm	2 cm	9
709163	16 G	15 cm	3 cm	9
7091635	16 G	15 cm	3,5 cm	9
709164	16 G	15 cm	4 cm	9
7091645	16 G	15 cm	4,5 cm	9
709362	16 G	20 cm	2 cm	9
709363	16 G	20 cm	3 cm	9
7093635	16 G	20 cm	3,5 cm	9
709364	16 G	20 cm	4 cm	9
7093645	16 G	20 cm	4,5 cm	9

FUSION 8 PLUS NEEDLES WITH SIDE INJECTION PORT AND THERMOCOUPLE:

CAT. NO.	CANNULA DIAMETER	LENGTH	ARRAYS DIAMETER	NUMBER OF ARRAYS
708162	16 G	15 cm	2 cm	8+1 Thermocouple
708163	16 G	15 cm	3 cm	8+1 Thermocouple
7081635	16 G	15 cm	3,5 cm	8+1 Thermocouple
708164	16 G	15 cm	4 cm	8+1 Thermocouple
701645	16 G	15 cm	4,5 cm	8+1 Thermocouple
708362	16 G	20 cm	2 cm	8+1 Thermocouple
708363	16 G	20 cm	3 cm	8+1 Thermocouple
7083635	16 G	20 cm	3,5 cm	8+1 Thermocouple
708364	16 G	20 cm	4 cm	8+1 Thermocouple
7083645	16 G	20 cm	4,5 cm	8+1 Thermocouple

SINGLE TIP NEEDLES

CAT. NO.	CANNULA DIAMETER	LENGTH	ACTIVE TIP
7013705	17 G	20 cm	0,5 cm
701371	17 G	20 cm	1 cm

RF GENERATOR

CAT. NO.	DESCRIPTION
MSI S-1500	150W RF GENERATOR, DUAL MODE CONTROL

COAXIAL SETS

CAT. NO.	DESCRIPTION
758152	COAXIAL INTRODUCER SET 12,5 CM
758157	COAXIAL INTRODUCER SET 17,5 CM

