

UNBLOCK

Fallopian Tube Set



CE 1936

HSG and Fallopian Tube Set

- 6 Fr catheter with 15° angled tip and transparent straightener
- Nitinol guidewire, with radiopaque polyurethane surface and hydrophilic coating



FEATURES



Radiopaque, double lumen catheter, with white and green color-coded connectors and 1-way FLL stopcocks, respectively for the balloon inflation and the contrast media infusion



The catheter tip has a soft balloon, which helps the anchoring within the uterus and prevents the reflux of the contrast media



The MLL syringe, connected to the stopcock, allows to keep the balloon volume and helps to deflate it quickly at the end of the procedure.

Ref. Number	CATHETER	GUIDEWIRE
MB855015	6 Fr, 40 cm (working length 27 cm)	0.035", 70 cm

