

RELIEF

Pleural and
Peritoneal Set

- 1 Safety valve
- 2 Polyester cuff
- 3 15.5 Fr silicone catheter
- 4 Fenestrated lenght 26 cm



RELIEF Pleural and Peritoneal Drainage Set

The Relief System is the ideal approach to the home management of recurrent pleural effusion and ascites and *is an effective alternative treatment to the procedures currently in use.*



FEATURES

- Safe and effective
- Easy to use
- Female Luer lock valve for universal connection
- Improved quality of life and reduced hospitalization times
- Active vacuum system drainage
- No complications
- Immediate improvement in respiratory performance status

Luer Lock valve

Prevents the accidental passage of air or fluid through the catheter and can be connected to any collection device through a MLL

9mm Dacron cuff

Minimizes the infection risks and avoids accidental displacement of the catheter.

Silicone catheter

15.5 Fr Soft and flexible catheter, long-term material with reduced kinking risk.

Barium Sulfate Line

Runs the entire length of the catheter, for an excellent X-ray visualization

25,5 cm Fenestrated portion

30 round holes ensure maximum drainage

Cat. No. MB 856016 The Relief Pleural and Peritoneal Drainage Kit contains:

Description	Quantity
Relief pleural catheter 15.5 Fr, 65 cm with FLL valve	1
Peel-away introducer 16 Fr, 17 cm	1
Tunneler 15,5 cm	1
Scalpel # 11	1
Stainless steel guidewire, FC, J3, .035", 45 cm	1
Needle 18 G, 9 cm	1
MLL Syringe 10 ml	1
Gauze pads	5
2000 ml graduated drainage bag with MLL connection	1
MLL/MLL Extension line, 100 cm	1
Drape 75 x 75 cm	1
Drape 100 x 100cm with adhesive hole diam. 20 cm	1
Foam Pad 10 x 10 cm	1
FLL bag adapter	1
Valve cap	1

Accessories:

- Cat. No. MB 807820 Drainage Vacuum Bottle 600 ml with MLL extension line
 Cat. No. MB 856017 Drainage Vacuum Bottle 600 ml Complete Set
 Cat. No. MB 807842 Drainage bag 2000 ml with MLL extension line (box of 25)
 Cat. No. MB 899002 Valve cap (box of 25)

